

JAY XTREME ACTIVE

Increased skin protection for active users



HIGH SKIN PROTECTION - LIGHTWEIGHT, LOW PROFILE

The JAY Xtreme Active cushion is a lightweight and low profile cushion for active performance and easy transfers. Designed to provide both high pressure relief and stability to active users, without compromising on weight or ride height. The ideal solution for active users.

The **JAY fluid pad** maximises skin protection and reduces shear.
The **dual cover system** promotes air flow and dissipates heat and moisture effectively.
Thigh positioning can be achieved through retro-fittable **leg support components**.





THE MODULARITY OF THE XTREME ACTIVE CUSHION

The JAY Xtreme Active cushion offers two different fluid options.

Standard fluid pad: Developed from the JAY Xtreme base and fluid pad, the standard JAY fluid insert offers moderate to high levels of pressure redistribution for most users, and goes with the standard base.

Large fluid pad: Developed from the JAY Active base and fluid pad, the larger option is designed for users with advanced muscle atrophy who do not want to go to a heavier high profile cushion but require more fluid under their bony areas. The large fluid pad has a higher fluid amount for greater immersion and pressure redistribution around the prominent bony areas of the pelvis and hips. It goes with a base with larger pelvic load area.



→ **Standard fluid pad**
With standard base.



→ **Large fluid pad**
In case of advanced muscle atrophy. It comes with a base with a larger pelvic load area.



LEG SUPPORT COMPONENTS

The Xtreme Active offers optional field-adjustable leg support components. The 2.5cm high adductor and abductor components can be easily placed on the inner cover using velcro. The thighs can be properly positioned to encourage orthopaedic alignment and increase seating tolerance.

→ **Leg support kit**
Contains two adductors
and one abductor.



Outer cover

Inner cover

Lightweight foam base



→ Microclimatic cover



→ Stretch cover



→ Incontinence cover

DUAL-COVER TECHNOLOGY

- The JAY Xtreme Active cushion features a dual-cover system. This innovative cover technology of an inner and outer cover provides an excellent microclimate and effective incontinence management.
- The inner cover material is water-resistant and easy to wipe off and clean, whilst its Aqua-guard zipper and anti-wicking thread ensure that the foam base stays dry. At the same time the microclimatic outer cover utilises the 3DX™ fabric to dissipate heat and moisture to maintain skin integrity.
- The JAY Xtreme Active is also available with an incontinence cover or outer stretch cover for greater immersion and envelopment into the fluid. All cushions come with a useful accessory bag stitched to the cover.
- Each inner and outer cover fits both fluid pad sizes and allows the leg supports to be fitted without creating surface tension.
- All covers are machine washable at 60°C and quick to dry.





TECHNICAL DATA

Cushion width	34cm – 50cm (in 2cm increments)
Cushion depth	34cm – 50cm (in 2cm increments)
Cushion height	5cm
User weight capacity	150kg
Cushion weight	1.7kg (40 x 40 with standard fluid pad); 2.2 kg (40 x 40 with large fluid pad)
Inserts	Standard & large JAY Flow Fluid™
Outer covers	Microclimatic, stretch and incontinent
Leg support kit	Contains 2 adductors & 1 abductor; each 2.5cm high (optional)

Find out more today!

Visit www.SunriseMedical.co.uk for further information on product options, or contact your local dealer or representative.



Get SOCIAL with QUICKIE and JAY

Keep up-to-date with the latest news on Quickie and JAY products through our Facebook page and YouTube channel:



www.facebook.com/QuickieWheelchairsUK

www.youtube.com/quickiewheelchairseu



Sunrise Medical Limited

Thorns Road, Brierley Hill, West Midlands, DY5 2LD

Tel.: +44 (0) 845 605 66 88 Fax: +44 (0) 845 605 66 89

E-mail: enquiries@sunmed.co.uk

www.SunriseMedical.co.uk