

# XENON<sup>2</sup>

  
**QUICKIE**<sup>®</sup>  
Live without limits.



**SUPERIOR FOLDING  
AND  
DRIVING  
PERFORMANCE**

# XENON<sup>2</sup> – THE STYLISH, HIGH-PERFORMANCE SERIES

---

Whether you're looking for the very best in driving dynamics or the most compact folding dimensions, there's a Xenon<sup>2</sup> just for you.

The New Xenon<sup>2</sup> is available in three different versions:

- The streamlined Xenon<sup>2</sup> fixed front
- The reinforced Xenon<sup>2</sup> hybrid fixed front
- The portable Xenon<sup>2</sup> swing away

Whichever one of the three feature packed models you choose, the sophisticated elegance of the Xenon<sup>2</sup> will always make you the centre of attention.





## XENON<sup>2</sup> FF THE LIGHTEST

The Lightest folding wheelchair in the world. Open fixed front frame with a cool, clean and streamlined design.



## XENON<sup>2</sup> H (HYBRID) THE STRONGEST

The hybrid is reinforced and has a fixed front frame, no rigidity is compromised and it also has a maximum user weight of 140 kg.



## XENON<sup>2</sup> SA THE ALLROUNDER

The Xenon<sup>2</sup> SA has the most compact folded dimensions for easy portability by car. The frame is reinforced and has a swing away front, for easier standing transfers in and out of your wheelchair.







# THE LIGHTEST FOLDING CHAIR IN THE WORLD...

---

How do you achieve a folding wheelchair that weighs as little as 8.8 kg? By using premium frame materials like aerospace-grade aluminium and an innovative process that's unique only to Quickie – ShapeLoc Technology.

By heat-treating the aluminium, we're able to increase its strength, which allows for the frame's tube walls to be thinner. Combine this with an oval design and hydroformed ShapeLoc Technology, this gives us a lighter and stronger wheelchair than anyone else. Just feel it...

# ...WITH UNPARALLELED FOLDING PERFORMANCE

---

At the core of the Xenon<sup>2</sup> is a new and unique cross-brace that fits so neatly under the seat it's barely noticeable. This creates a very minimalistic, open-frame feel that you'd usually associate with a rigid wheelchair.

It doesn't just end there! The new axle stem has been completely re-designed and combined with the new cross-brace, provides the stiff driving performance and feel of a rigid frame.

## Easy to use folding mechanism

The new aluminium latches on the cross-brace will also keep your wheels in parallel when folding, so your Xenon<sup>2</sup> remains stable no matter what the situation.






# STYLE IN MOTION

---

Personalise your Xenon<sup>2</sup> and stand out from the crowd with a versatile palette of 32 frame colours. There's also 5 distinctive anodised styles to choose from for the castor fork and castor wheel. You can even accentuate your colour scheme further by matching or contrasting the upholstery trim colour of the EXO backrest upholstery.

→ Personalise your Xenon<sup>2</sup>  
With a selection of 32  
eye-catching colours







→ Add some colour to your back upholstery  
This very comfortable and stylish Exo back upholstery can be chosen with 5 different colour trims.



→ Choose your own castor and fork style  
The anodised front castors and the forks are available in silver, blue, red, orange or black for a unique look.





# ADJUSTABLE TO SUIT YOU

---

Exactly as it should be – designed to be fine-tuned and adaptable to your changing needs.

Do you want to change your chair driving characteristic from passive to more active, or the opposite way? Just change the centre of gravity setting by simple repositioning of the axle stem: and don't forget to adjust the back angle for more ergonomic posture.

- 1 **Centre of Gravity**  
Enables a driving characteristic from very active to passive.
- 2 **Backrest Angle**  
The right back angle adjustment is important for your ergonomic posture.
- 3 **Rear Seat Height**  
Exactly as it should be - adjust the height to your personal need.





**XENON² FF**  
**THE LIGHTEST**



**XENON² H (HYBRID)**  
**THE STRONGEST**



**XENON² SA**  
**THE ALL ROUNDER**

# **XENON²**

Ultra lightweight wheelchairs  
with outstanding folding and  
driving performance





Options and features of the Xenon<sup>2</sup> series – Configure your wheelchair to your individual needs. For further information on product options, please check the Xenon<sup>2</sup> and Xenon<sup>2</sup> SA order form.



**// FOLD-DOWN BACKREST**

This half-folding backrest folds to the front and achieves a small package size, making it easier to transport.



**// SURGE® LT HANDRIMS**

Its Gripton rubber strip on top of the ergonomic oval shaped rim offers maximum traction.



**// ANGLE ADJUSTABLE BACK**

Stylish designer backrest solution integrated into the back of the frame. Backrest angle from 75° to 103° in 4° increments.



**// DESK SIDEGUARD**

This flip-up and removable sideguard can be ordered with fixed or height adjustable armpads.



**// PADDED TUBULAR ARMREST**

These armrests give you the support when needed and do not hinder during propulsion, they can also swing away.



**// CARBON SIDEGUARD**

Very strong and lightweight sideguard with a unique look.



**// AUTO-FOLDING FOOTBOARD**

This footboard is angle and depth adjustable and available with carbon or composite. Only available for fixed front frame versions.



**// TIE-DOWN BRACKETS FOR TRANSPORTATION**

With this option the wheelchair is approved for transportation in a vehicle (crash tested according to ISO 7176-19).



# TECHNICAL DATA

	Xenon <sup>2</sup> FF	Xenon <sup>2</sup> H (Hybrid)	Xenon <sup>2</sup> SA
Max. user weight:	110 kg (17 stone)	140 kg (22 stone)	125 kg (20 stone)
Seat width:	320 – 460 mm	320 – 500 mm	320 – 500 mm
Seat depth:	340 – 500 mm	340 – 500 mm	340 – 500 mm
Seat height front:	430 – 550 mm	430 – 550 mm	430 – 560 mm
Seat height rear:	370 – 500 mm	370 – 500 mm	370 – 500 mm
Backrest height:	250 – 475 mm	250 – 475 mm	250 – 475 mm
Backrest angle:	-15° to +13° increments	In 4° increments	In 4° increments
Centre of gravity (COG):	40 – 130 mm	40 – 130 mm	40 – 130 mm
Camber:	0° / 2° / 4°	0° / 2° / 4°	0° / 2° / 4°
Frame Version:	Open frame with fixed front	Reinforced frame with fixed front	Reinforced frame with swing-away front
Chair weight:	From 8.8 kg	From 9.5 kg	From 10.3 kg
Frame colours:	32	32	32
Crash tested:	Yes	Yes	Yes

For further information on the full specification, options and accessories, please check the order form on our website.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk) or alternatively it is available on request in large text.



[www.facebook.com/quickiewheelchairs](http://www.facebook.com/quickiewheelchairs)

**Sunrise Medical Limited**

Thorns Road, Brierley Hill, West Midlands, DY5 2LD

Tel: +44 (0) 845 605 66 88 Fax: +44(0) 845 605 55 89

E-Mail: [enquiries@sunmed.co.uk](mailto:enquiries@sunmed.co.uk)

[www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk)