

Q300 M Mini Teens

The ultra-slim powerchair
you won't grow out of.



THE NARROWEST CHAIR FOR TEENS ON THE GO (AND ON THE GROW).

It's a fact of life. Teenagers grow. Fast. Faster than their clothes and faster than their powerchair. That's why we created the **Q300 M MINI Teens** - a powerchair to keep up with growing teens!

With an adjustable seating system that grows by height and width, getting older no longer means a big adjustment to a new chair. Instead, make lots of small ones and your Q300 Mini Teens will fit you at 18 just as it did at 13!

But of course it wouldn't be mini without a ridiculously narrow, super-manouverable **520 mm** base to get you out of tight situations. Say goodbye to Dad-Taxi too. From school, cinema to hanging out with friends, it's packed with all the big outdoor performance you need for an independent life (and with none of the size).

Plus, with **C-Me** you can drive with your seat elevated at walking speed, so you can keep up with your friends (both at eye-level and pace) too!



Scan the QR code with your phone's camera to watch video.



The Q3 narrow or tight





Q300 M Mini Teens is so ridiculously
small, squeezing between parked cars
and doorways is a breeze!



520

MM
WIDTH

TRUE MID-WHEEL-DRIVE TECHNOLOGY
The Q300 M Mini Teens can provide the tightest
turning circle with increased traction and stability
thanks to its mid-wheel drive technology.



1120 mm
COMPACT TURNING CIRCLE

COMPACT SPACE? HELLO MID-WHEEL DRIVE BASE!

From cluttered classrooms to hectic hallways, the Q300 M MINI Teens is so ridiculously compact and manoeuvrable, it'll easily fit into your daily life (not the other way around).

410 mm

SEAT-TO-FLOOR HEIGHT

A real low-rider, the Q300 M Mini Teens effortlessly fits under those classroom tables with an ultra low seat-to-floor height starting from **410 mm** (with fixed tilt) and **430 mm** (with powered lift/tilt).



CAPSULATED JOYSTICK

Joystick not in the right position? This (optional) mounting option has 3 locking positions - including **outwards** (for getting closer to tables), **inwards** (if you have limited hand mobility) and the **standard** driving position.

1120 mm

BASE LENGTH

Ultra-short in length, the Q300 M Mini Teens conveniently fits in-and-out of the smallest lifts and tight spaces with ease. It's also light enough (from **130 kg**) to be compatible with most floor lifts too.

GET CLOSER WITH THE CENTRE MOUNT LEGREST

Swap traditional legrests at the front of your powerchair with the (optional) centre-mount and not only can you get closer to objects (like under tables), but you'll get the tightest turning circle too!



BIG POWERCHAIR PERFORMANCE (MINUS THE BIG)

Now this is one compact powerchair with BIG outdoor performance. Narrow and nimble in a crowded city, and packing a ton of outdoor features, it's streets ahead of the competition.

70 mm KERB CLIMBING

No drop-kerb in sight? The Q300 M MINI Teens climb kerbs of up to **70 mm**, so there's no need to take the long-way around.



36 km
RANGE

Powered by **56 Ah batteries**, you'll get everywhere you want to go without having to rush back with a range of up to **36km** (22 miles).

can



GET TO GRIPS WITH GREATER TRACTION

Get the best traction, stability and obstacle tackling with outdoor **14"** drive wheels. They'll also keep your Q300 M MINI slim at just **570 mm**.



PUSH PERFORMANCE WITH 4-POLE MOTORS

These motors provide a punchier response to driving commands and also a greater ability to tackle troublesome terrain. Longer-lasting, highly efficient with improved performance, our motor technology can't be easily duplicated anywhere else.

ANTI-PITCH TECHNOLOGY

Traditionally when a powerchair comes to a halt on a slope, the stopping motion will jolt you forwards. We've solved that problem with Anti-Pitch technology, so you're kept safe and stable in your seat.



ALL-WHEEL INDEPENDENT SUSPENSION

With our **PATENTED** all-wheel independent suspension, all six wheels remain on the ground whilst climbing and transitioning obstacles. As any shocks or bumps are absorbed by the wheel's suspension, you're provided a soft and smooth ride.

SEATING THAT GROWS WITH YOU.

Growing up doesn't mean growing out of your powerchair. We've designed the Q300 M MINI Teens seat so you never have to fear those growth spurts. Adjustable by width (280-400 mm) and depth (360-500 mm), it's ideal for those requiring a narrow seat width, but greater adjustments to depth over time.



UNIVERSAL HEADREST
Compatible with **SEDEO** and **WHITMYER** headrests, this headrest mount will help to address a range of head-support needs (from basic to multiple-point control).

Z-FINITY FOOTRESTS

Adjustable between **180 mm** above the seat pan to **480 mm** below, these footrests provide support for growing teens. Its pivot arm allows for infinite adjustment of the footrest angle and placement (for ideal positioning of the feet, ankles and knees).



SEDEO LITE DEEP CONTOUR BACKREST

Deeply-contoured for comfort, this backrest has tension-adjustable straps that can be easily adjusted to follow and support a variety of postural shapes. The deep laterals are also integrated inside the backrest, so not only will they provide constant trunk stability, they also won't move or slip over time either.



DO YOU HAVE SPECIFIC CLINICAL SEATING REQUIREMENTS?

The Q300 M MINI Teens seat is also compatible with the **J3® Back**. Suitable for moderate-to-high positioning, the J3® Back offers the widest range of sizes, shapes and options to precisely fit and tailor your seat just to you.

POWER-UP YOUR POWERCHAIR!

The Q300 M MINI Teens doesn't just have BIG chair performance, it's also loaded with BIG features to meet your clinical and lifestyle needs.

30°

TILT-IN-SPACE

Adjust your seated position throughout the day with 30° of powered tilt. Helping with pressure relief by redistributing your weight from the cushion to the backrest, you can maintain a comfortable seated position with reduced muscular effort.



R-NET ADVANCED JOYSTICK

Ramp up the on-board technology with this (optional) large, easy-to-read HD display and simple to control paddle switches. Take control of your world with integrated Bluetooth and i-device functionality and control your TV, phone and other devices all from your Q300 M Mini Teens.

POWERED ELEVATING & ARTICULATING LEGRESTS

By raising and lowering your legs, these legrests help to reduce swelling and improve circulation. By adjusting in length, your seated position is kept constant, preventing unwanted pelvic movement.

POWERED LIFT & TILT

Reach new heights with **300 mm** of seat lift and get to everyday objects that would otherwise be out of grasp. You can also transfer easier thanks to being in a higher position to stand from. Plus, with 30° of tilt, you can sit back and enjoy greater relaxation and pressure distribution too.





ELEVATE AND DRIVE WITH C-ME

Say goodbye to waiting around or driving at reduced speed when elevated. With C-Me, you can raise up to 300 mm and drive up to 4.8 kph all at the same time! Now you can keep up with your friends at their pace and be at eye-level all with a touch of a button (with no need to stop).



RAISING COMFORT WHEN DRIVING ELEVATED

Unlike other powerchairs, our suspension doesn't lock-out when using C-Me. So you get the same smooth and comfortable ride over uneven terrain, even when elevated! Plus, anti-pitch technology will keep you safe and stable on inclines too.





STOP SWITCH FOR ADDED SAFETY

There's safety in stopping. Our added remote stop switch allows an attendant to stop a powered wheelchair via an RF remote from up to 30 meters away. Operating distance dependent on the environment.



ELEVATE WHILST DRIVING WITH ASSIGNABLE BUTTONS

Fed up of complex joystick movements or having to stop to use your seat lift? With our **PATENTED** assignable buttons, you can change any button (or assign a second, long-press function) to activate C-ME (even when driving).




BRIGHT YOUNG THING.

You'll be sitting pretty with 12 fresh shroud designs!
Co-ordinate or contrast them with 7 accent colours
for the wheel and castors to get the perfect mix.



SHROUD:



ACCENT:



Want something unique? Built-4-Me can help you create a truly bespoke Q300 M MINI Teens that's really out of this world!





ONE-OF-A-KIND UPHOLSTERY COLOURS

Now this is upholstery that's really cut from a different type of cloth! Our Built-4-Me team can personalise your seat with a variety of unique upholstery colours and really turn heads!



PERSONALISED SHROUD DESIGNS

Prefer a more individual feel to your powerchair shroud? Almost any type of design can be printed and foil-wrapped to your shroud by our Built-4-Me team.

TECHNICAL DATA

Seat width	280 - 400 mm (11 - 16")
Seat depth	360 - 500 mm (14 - 20")
Seat height	410 - 470 mm (16 - 19")
Back height	385 - 510 mm (15 - 20")
Lower leg length	180 - 510 mm (7 - 20")
Powered tilt	0 - 50° / 0 - 30° / -5 - 25°
Total width	From 520 mm
Total length (with centre-mount)	1060 mm
Turning circle	1120 mm
Maximum safe slope	6°
Maximum kerb climbing	50 mm (with 12" drive wheels) 70 mm (with 14" drive wheels)
Maximum range (ISO 7176-4)	Up to 25 km (with 41 Ah batteries) Up to 36 km (with 56 Ah batteries)
Product weight (with batteries)	From 103 kg
Maximum user weight	120 kg
Drive wheels	12" / 14"
Maximum speed	6 kph / 8 kph / 10 kph
Battery size	41 Ah / 50 Ah / 56 Ah
Electronics	VR2, R-NET
Crash tested	Yes - approved for forward facing transport in a vehicle (ISO 7176-19). Attachment: 4 chair fixing points with 2 double straps rear and 2 single straps front. No anchored belt restraints allowed. The ease of access and manoeuvrability in vehicles depends on the wheelchair size.
Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane).
Intended use and environment	For people with limited mobility, indoor and outdoor use. Wheelchair type class B



Scan the QR code with your phone's camera and book your free test-drive today.

Now that you've read all about the QUICKIE Q300 M Mini Teens, why not arrange a free demonstration and experience it for yourself?

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines.

In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit www.sunrisemedical.co.uk

If you are visually impaired, this document can be viewed in PDF format at www.sunrisemedical.co.uk



Sunrise Medical Limited

Thorns Road, Brierley Hill, West Midlands, DY5 2LD

Tel.: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

E-Mail: enquires@sunmed.co.uk

www.sunrisemedical.co.uk

