

# Q300<sup>®</sup> M Mini Kids

Imagine the Possibilities



# ZIPPIE Q300 M MINI KIDS

Whether it's navigating round the home, school or exploring outdoors, the Q300 M Mini Kids has just as much energy for adventure as your child!

The Q300 M Mini Kids with its 520 mm **ultra narrow base and small turning circle** of 1120 mm is made for easily and intuitively moving around in tight spaces or discovering out-of-the-way corners outdoors. And when adventures take them further than you thought possible, the **super-strong 4-pole motors** (optional) will power through. At the same time, the **all-wheel suspension** ensures a smooth and safe ride when things get a bit wild.

Together with the **adjustable seating system** that makes sure that the kids don't grow out, there are almost endless styling possibilities according to individual taste. With helpful powered options the Q300 M Mini Kids is designed for your child's independence.



Scan the QR code with your phone's camera to get more information.



## Additional Safety Feature

If the little ones are too fast and far away, you can use the Remote Switch Stop to stop the wheelchair from 30 m away (depending on environment).

# IMAGINE THE PLACES YOU'LL GO

Let your child discover the world the way he or she wants to and let it experience the most amazing adventures!

The mid-wheel drive technology allows most intuitive and easy driving for your child. With a drive wheel that's positioned in the middle of the base, your kid is able to turn on the spot for great indoor manoeuvrability. Plus, as the child is seated directly over the drive wheel, the positioning of it's weight increases traction to provide exceptional outdoor performance too.

Looking for the most compact base for indoor manoeuvrability? Choose the 12" drive wheels, if you prefer better outdoor performance? Go for 14"!



Mid-wheel drive technology



### All-wheel independent suspension

Our PATENTED all-wheel independent suspension ensures all six wheels remain on the ground whilst climbing and transitioning obstacles. Any shocks or bumps are absorbed by the wheel's suspension, providing a soft and smooth ride.

### Anti-Pitch technology

Traditionally when a powerchair comes to a halt on a slope, the stopping motion will jolt you forwards. That's a problem we've eliminated with Anti-Pitch technology - so your child kept safe and stable in it's seat.

### Kerb climbing

The Q300 M Mini Kids with its 300 mm (12") drive wheels can handle kerbs of up to 50 mm so there's no need to take the long way round. Upgrade to 355 mm (14") drive wheels and tackle kerbs of up to 70 mm.

# A SEATING SYSTEM THAT GROWS WITH YOUR CHILD

**Children grow, and sometimes very quickly.**

With the Q300 M Mini Kids, there's no need to worry about to changing to a new chair after your kids have just got used to the one wheelchair since the seating system has been specifically designed to grow with your child.

The **adjustable seating system** can grow without additional parts in width (280 - 400 mm) and depth (260 - 400 mm) and makes it fit properly. This leads to the best positioning and good stability and makes your child feels comfortable throughout the day, both while playing as well as during the necessary rest periods.

This makes the Q300 M Mini Kids the companion for all small and big adventures in life.



Scan the QR code with your phone's camera to get more information.



## **Z-Finity Footrest**

Infinitely adjustable (between 160 mm above the seat pan to 480 mm below) along its rotational arc, multiple angles and placement options for positioning the feet, ankles, and knees can be achieved.



### Powered Lift and Tilt

Reach new heights with 300 mm of seat lift and get to every-day objects that would otherwise be out of grasp. Your child can also transfer easier thanks to being in a higher position to stand from. Plus, with up to 50° of tilt, your child can sit back and enjoy greater relaxation and pressure distribution too.

### Universal Headrest

The headrest mount on Q300 M Mini Kids is compatible with SEDEO and WHITMYER headrests, which help to address a range of needs from basic support to multiple-point control.



### SEDEO LITE Backrest

We've designed this backrest just for kids. Deeply contoured for comfort, the integrated laterals provide constant trunk stability. As they're inside the backrest, the laterals won't move or slip over time - great for kids on the move! And with eight tension-adjustable straps, a variety of postural shapes can be supported too.



### Capsulated Joystick

This joystick mount has three locking positions. It can be mounted mid-line for children who have limited hand mobility and can also swing inwards/outwards to help getting closer to tables.



# ELEVATE AND DRIVE WITH C-ME

## C-Me is all about seeing your child, and just your child!

Drive up to 4.8 km/h while elevated up to 300 mm without needing to stop the seat position. Your child can drive shoulder to shoulder with family or friends.

Driving at walking speed in full lift while being together with the classmates or enjoying a stroll with the family - C-Me takes away frustrating speed restrictions allowing your child being in the middle of just being there.

Let your child discover the world the way he or she wants to, without being dependent on you or others.



## Ride in comfort

Our power wheelchairs are designed with style for long-lasting performance. Unlike some of our competitors, our suspension doesn't lock out when you drive with C-Me, giving you a smooth and comfortable ride from start to finish.



Scan the QR code with your phone's camera to get more information.

# BRIGHT YOUNG THING



## You want to sit in style?

You will be with 12 exciting shroud designs patterns!  
Coordinate or contrast them with an additional seven accent colours for the wheel and casters to get the perfect mix.



### SHROUD:



### ACCENT:



## Looking for something unique?

A person's individuality is as unique as their fingerprint. That is why we are also happy to design wheelchairs according to your kid's personal wishes - like a personally designed foil wrapping. Considerably more detailed and versatile, almost any type of design can be printed onto foil and affixed to your wheelchair's frame, sideguards and even footrest.

## The only limitation is your imagination!

Individual wheelchair frame design foil



But there's more. Our BUILT-4-ME department has created thousands of bespoke powerchairs. When you want to break the mould, we'll help you create a truly one-of-a-kind powerchair.

Find out more at:

[www.sunrisemedical.eu/built4me](http://www.sunrisemedical.eu/built4me)

# TECHNICAL DATA



Scan the QR code with your phone's camera to get more information.

## Q300 M Mini Kids

<b>Seating System</b>	Kids seat	<b>Lower leg length</b>	180 - 460 mm (7" - 18")
<b>Powered tilt</b>	0° - 50° / 0° - 30° / -5° - 25°	<b>Seat height</b>	410 - 490 mm (16" - 19")
<b>Seat width</b>	280 - 400 mm (11" - 16")	<b>Max user weight</b>	75 kg
<b>Seat depth</b>	260 - 400 mm (10" - 16")	<b>C-Me</b>	Yes
<b>Back height</b>	380 - 485 mm (15" - 20")	<b>Maximum range (ISO 7176-4)</b>	Up to 25 km (with 41 Ah batteries) Up to 36 km (with 56 Ah batteries)
<b>Total width base</b>	From 520 mm	<b>Product weight (with batteries)</b>	From 103 kg
<b>Total length (with centre-mount) base</b>	1060 mm (41.7")	<b>Drive wheels</b>	12" / 13" / 14"
<b>Turning circle</b>	1120 mm (44.1")	<b>Maximum speed</b>	6 km/h / 8 km/h / 10 km/h
<b>Maximum safe slope</b>	6°	<b>Battery size</b>	41 Ah / 50 Ah / 56 Ah
<b>Maximum kerb climbing</b>	50 mm (with 12" drive wheels) 60 mm (with 13" drive wheels) 70 mm (with 14" drive wheels)	<b>Transport without user</b>	Fix the wheelchair safe in the vehicle (e.g. car or plane)
<b>Electronics</b>	VR2, R-NET	<b>Intended use and environment</b>	For people with limited mobility, indoor and outdoor use. Wheelchair type class B



### Sunrise Medical GmbH

Kahlbachring 2-4  
69254 Malsch/HD  
Tel.: +49 (0) 72 53 / 9 80-0  
Fax: +49 (0) 72 53 / 9 80-222  
E-Mail: [international@sunrisemedical.de](mailto:international@sunrisemedical.de)  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.eu](http://www.SunriseMedical.eu).

